

AériaTMMC

Orthokeratology lens
with Myopia Control

Intended use

Ultra permeable Orthokeratology lens:

- Corneal reshaping
- Myopia and astigmatism
 - Myopia up to - 10.00 D
 - Astigmatism up to 4.00 D
- Myopia Control

Designs

You can choose between to 2 designs:

- **AériaTMMC** : Myopia and corneal astigmatism up to 1.00 D
- **AériaTMMC_{Toric}** : Myopia and corneal astigmatism above 1.00 D


Technical specifications

Material	Dk 200 Dissociate Colors RE-LE
Diameter Ø _T	From 10.50 to 12.00 mm by step of 0.10 mm
Base Curve	From 7.50 to 11.00 mm by step of 0.10 mm
Reverse	From 10 to 200 µm by step of 10 µm
Toricity	From 20 to 100 µm by step of 10 µm
Peripheral alignment	
K reading	From 7.40 to 8.80 mm by step of 0.10 mm
e (corneal eccentricity)	From 0.3 to 0.8 by step of 0.1
Power P	From -2.00 to + 2.00 D by step of 0.25 D

Wear and Care

Wear	Annual renewal, overnight wear
Care	Ote Clean (Cleaner) and RGP Multipurpose solution
Handling	Insertion: instill 2 drops of moisturising solution in the lens at the insertion. Removal: like an RGP lens, with a suction cup

WWW.LABORATOIRE-LCS.COM



AériaTMMC

Orthokeratology lens
with Myopia Control


Calculate the first lens

To calculate the Aéria MC parameters, thank to send us the information below (export@laboratoire-lcs.com):

- Corneal diameter
- Spectacle refraction
- Pupil diameter
- Keratometry
- Corneal Eccentricity] or the corneal topography


Fitting assessment

1. Centration




2. Fluo Control

Two peripheral rings of reservoir are specifically placed in peripheral zone to help you to control the well alignment to control the well alignment




Fluo under
the internal reservoir




Open peripheral alignment
↗ eccentricity

No fluo under the both
peripheral reservoirs



Optimal
Alignment

Fluo under
the external reservoir



Close peripheral alignment
↘ eccentricity

export@laboratoire-lcs.com



WWW.LABORATOIRE-LCS.COM

LCS © ALL RIGHTS RESERVED 2020 - 06/20